



Product Datasheet

GKN1 antibody (orb537823)


Description	Rabbit polyclonal antibody to GKN1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GKN1(NP_062563.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22kDa
Uniprot ID	Q9NS71
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human stomach tissue.



Immunohistochemical staining of human placenta tissue.



Western blot analysis of extracts of mouse stomach.