



Product Datasheet

NGAL antibody (orb537844)


Description	Rabbit Polyclonal antibody to NGAL
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human Lipocalin-2/NGAL(NP_005555.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22kDa
Uniprot ID	P80188
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse kidney...



Immunofluorescence analysis of rat kidney...



Western blot analysis of extracts of mouse...