


## Product Datasheet

### SUN1 antibody (orb537854)

|                            |                                                   |
|----------------------------|---------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SUN1                |
| <b>Reactivity</b>          | Human, Mouse, Rat                                 |
| <b>Conjugation</b>         | Unconjugated                                      |
| <b>Tested Applications</b> | WB                                                |
| <b>Immunogen</b>           | Recombinant fusion protein of human SUN1(O94901). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.  |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.       |
| <b>Note</b>                | For research use only                             |
| <b>Isotype</b>             | IgG                                               |
| <b>Clonality</b>           | Polyclonal                                        |
| <b>MW</b>                  | 90kDa                                             |
| <b>Uniprot ID</b>          | <b>O94901</b>                                     |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                   |



Western blot analysis of extracts of mouse testes.