


## Product Datasheet

### SREBF1 antibody (orb537872)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SREBF1                      |
| <b>Species/Host</b>        | Rabbit                                                    |
| <b>Reactivity</b>          | Mouse                                                     |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Tested Applications</b> | IF, WB                                                    |
| <b>Immunogen</b>           | A synthetic peptide of human SREBF1 (NP_001005291.1).     |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.               |
| <b>Note</b>                | For research use only                                     |
| <b>Isotype</b>             | IgG                                                       |
| <b>Clonality</b>           | Polyclonal                                                |
| <b>Purity</b>              | Affinity purification                                     |
| <b>MW</b>                  | 125KDa                                                    |
| <b>Uniprot ID</b>          | <b>P36956</b>                                             |
| <b>Entrez</b>              | <b>6720</b>                                               |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |



Western blot analysis of extracts of mou...