

Product Datasheet

GAB1 (Phospho-Y659) antibody (orb537889)

Description Rabbit polyclonal antibody to GAB1 (Phospho-Y659)

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen A synthetic phosphorylated peptide around Y659 of human GAB1(NP_002030.2).

Preservatives PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 110kDa

Uniprot ID [Q13480](#)

Dilution Range WB: 1:500-1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of NIH...