


Product Datasheet

PDGFRb (Phospho-Y740) antibody (orb537895)

Description	Rabbit polyclonal antibody to PDGFRb (Phospho-Y740)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y740 of human PDGFRB (NP_002600.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	180kDa
Uniprot ID	P09619
Entrez	5159
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
extracts of
NIH...