



Product Datasheet

ACLY (Phospho-S455) antibody (orb537900)


Description	Rabbit polyclonal antibody to ACLY (Phospho-S455)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S455 of human ACLY (NP_001087.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	121kDa
Uniprot ID	P53396
Entrez	47
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse pancreatic islets



Immunohistochemical staining of mouse brain tissue



Immunohistochemical staining of rat pancreatic islets