


Product Datasheet

CDKN1B antibody (orb537907)

Description	Rabbit polyclonal antibody to CDKN1B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human CDKN1B/p27KIP1(NP_004055.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	P46527
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of PC-...