



Product Datasheet

APOBEC3G antibody (orb537917)


Description	Rabbit polyclonal antibody to APOBEC3G
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human APOBEC3G(NP_068594.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45KDa
Uniprot ID	Q9HC16
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cells...



Immunohistochemical staining of rat testis...



Western blot analysis of extracts of PC-3 cells...