



Product Datasheet

CD79A antibody (orb537921)

Description	Rabbit polyclonal antibody to CD79A
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human CD79a(NP_001774.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40-50kDa
Uniprot ID	P11912
Dilution Range	WB: 1:500-1:1000, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of Raj...