



Product Datasheet

HLA-DRB3 antibody (orb537932)


Description	Rabbit polyclonal antibody to HLA-DRB3
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human HLA-DRB3(NP_072049.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30kDa
Uniprot ID	P79483
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ap...



Immunohistochemical staining of human lu...



Immunohistochemical staining of tonsil t...