


Product Datasheet

BRAF (Phospho-S445) antibody (orb537939)

Description	Rabbit polyclonal antibody to BRAF (Phospho-S445)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S445 of human BRAF(NP_004324.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	100kDa
Uniprot ID	P15056
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Raj...