


Product Datasheet

PLCG2 (Phospho-Y759) antibody (orb537941)

Description	Rabbit polyclonal antibody to PLCG2 (Phospho-Y759)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y759 of human PLC gamma 2 (PLCG2)(NP_002652.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	150kDa
Uniprot ID	P16885
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Raj...