


Product Datasheet

KCNMA1 antibody (orb537972)

Description	Rabbit polyclonal antibody to KCNMA1
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human KCNMA1(NP_001258447.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	137kDa
Uniprot ID	Q12791
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Rat...