



Product Datasheet

KCNS2 antibody (orb537973)


Description	Rabbit polyclonal antibody to KCNS2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human KCNS2(NP_065748.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54kDa
Uniprot ID	Q9ULS6
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse eye...



Immunofluorescence analysis of mouse eye...



Immunofluorescence analysis of rat eye l...