

Product Datasheet

PLCB3 antibody (orb537981)

Description	Rabbit polyclonal antibody to PLCB3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human PLC beta 3 (PLCB3)(NP_001303243.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	160kDa
Uniprot ID	Q01970
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.