

Product Datasheet

STX5 antibody (orb537988)

Description	Rabbit polyclonal antibody to STX5
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human STX5(NP_003155.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q13190
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.