


Product Datasheet

LAMB3 antibody (orb538038)

Description	Rabbit polyclonal antibody to LAMB3
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human LAMB3(NP_000219.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130KDa
Uniprot ID	Q13751
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of rat...