



Product Datasheet

PNLIPRP2 antibody (orb538080)


Description	Rabbit polyclonal antibody to PNLIPRP2
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human PNLIPRP2(NP_005387.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54kDa
Uniprot ID	P54317
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse pancreas



Immunofluorescence analysis of rat pancreas



Western blot analysis of extracts of rat...