

Product Datasheet

Lysozyme-like protein 1 antibody (orb538108)

Description Rabbit Polyclonal antibody to Lysozyme-like protein 1

Reactivity Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Recombinant fusion protein of human LYZL1(NP_115906.3).

Preservatives PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 17kDa

Uniprot ID [Q6UWQ5](#)

Dilution Range WB: 1:500-1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of rat...