


Product Datasheet

PLEKHG4 antibody (orb538111)

Description	Rabbit polyclonal antibody to PLEKHG4
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PLEKHG4(NP_001123199.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	131kDa
Uniprot ID	Q58EX7
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of rat...