



Product Datasheet

PRDM14 antibody (orb538112)


Description	Rabbit polyclonal antibody to PRDM14
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IHC, IP, WB
Immunogen	Recombinant fusion protein of human PRDM14(NP_078780.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	Q9GZV8
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, IP: 1:20-1:100, ChIP: 1:20-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 cells...



Immunofluorescence analysis of HeLa cells...



Western blot analysis of extracts of rat...