


Product Datasheet

EphA2 (Phospho-Y594) antibody (orb538137)

Description	Rabbit polyclonal antibody to EphA2 (Phospho-Y594)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y594 of human EphA2 (NP_004422.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	108kDa
Uniprot ID	P29317
Entrez	1969
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of SH-594 cells.