


Product Datasheet

GRIA1 (Phospho-S845) antibody (orb538138)

Description	Rabbit polyclonal antibody to GRIA1 (Phospho-S845)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S845 of human GluR1/GRIA1
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	90kDa
Uniprot ID	P42261
Entrez	2890
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of SH-...