



Product Datasheet

TNFRSF1B antibody (orb538177)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to TNFRSF1B |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, WB |
| Immunogen | Recombinant fusion protein of human TNFR2/TNFR2/TNFRSF1B(NP_001057.1). |
| Preservatives | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 65kDa |
| Uniprot ID | P20333 |
| Dilution Range | WB: 1:500-1:2000, IF: 1:50-1:200 |
| Expiration Date | 12 months from date of receipt. |



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of THP...