



Product Datasheet

CREB3 antibody (orb538215)

Description	Rabbit polyclonal antibody to CREB3
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CREB3(NP_006359.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Uniprot ID	O43889
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue showing CREB3 expression in neurons.



Western blot analysis of extracts of U-87MG cells showing CREB3 expression.