

Product Datasheet

SLC25A23 antibody (Biotin) (orb53822)

Description	SLC25A23 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein SCaMC-3 protein (1-188AA)
Target	SLC25A23
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BV35
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.