


Product Datasheet

VEGFB antibody (orb538248)

Description	Rabbit polyclonal antibody to VEGFB
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human VEGFB (NP_003368.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	32KDa
Uniprot ID	P49765
Entrez	7423
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of U-8...