



Product Datasheet

CD35 antibody (orb538258)

Description	Rabbit Polyclonal antibody to CD35
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human CD35/CD35/CR1 (NP_000642.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	270kDa
Uniprot ID	P17927
Entrez	1378
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of THP-1 cel...



Western blot analysis of extracts of U-9...