



Product Datasheet

IL2 antibody (orb538263)


Description	Rabbit polyclonal antibody to IL2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human IL2(NP_000577.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15KDa/16KDa/18KDa
Uniprot ID	P60568
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lu...



Immunohistochemical staining of human lu...



Immunohistochemical staining of human to...