



Product Datasheet

ABCC3 antibody (orb538277)

Description	Rabbit polyclonal antibody to ABCC3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ABCC3 (NP_003777.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	169kDa
Uniprot ID	O15438
Entrez	8714
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human prostate tissue.



Western blot analysis of extracts of various cell lines.