


Product Datasheet


Acetyl-Histone H4-K16 antibody (orb538299)

Description	Rabbit polyclonal antibody to Acetyl-Histone H4-K16
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic acetylated peptide around K16 of human Histone H4(NP_001029249.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	11KDa
Uniprot ID	P62805
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, ChIP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells

...



Western blot analysis of extracts of var...