



Product Datasheet

ADAMTS13 antibody (orb538338)

Description	Rabbit polyclonal antibody to ADAMTS13
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human ADAMTS13(NP_620594.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	154kDa
Uniprot ID	Q76LX8
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Western blot analysis of extracts of var...