



Product Datasheet

FANCF antibody (orb53835)

Description	FANCF antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Fanconi anemia group F protein (295-341AA)
Target	FANCF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NPI8
Dilution Range	WB:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293(lane 1),...



Western blot analysis of HEK293 whole ce...