



Product Datasheet

AK4 antibody (orb538384)

Description	Rabbit polyclonal antibody to AK4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human AK4 (NP_037542.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	25kDa
Uniprot ID	P27144
Entrez	205
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Western blot analysis of extracts of var...