



Product Datasheet

AKT1S1 antibody (orb538400)


Description	Rabbit polyclonal antibody to AKT1S1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human AKT1S1/PRAS40(NP_001092103.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Uniprot ID	Q96B36
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Immunohistochemical staining of human st...



Immunohistochemical staining of human ma...