


Product Datasheet

ANT1/ANT2/ANT3/ANT4 antibody (orb538433)

Description	Rabbit polyclonal antibody to ANT1/ANT2/ANT3/ANT4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human ANT1(NP_001142.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30kDa
Uniprot ID	P12235 , Q9H0C2 , P05141 , P12236
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...