



Product Datasheet

APLP1 antibody (orb538449)


Description	Rabbit polyclonal antibody to APLP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human APLP1(NP_005157.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100kDa
Uniprot ID	P51693
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue showing APLP1 expression in neurons.



Immunohistochemical staining of human cerebellar tissue showing APLP1 expression in neurons.



Western blot analysis of extracts of various cell lines showing APLP1 protein expression.