



Product Datasheet

AQP2 antibody (orb538465)


Description	Rabbit polyclonal antibody to AQP2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human AQP2(NP_000477.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26KDa
Uniprot ID	P41181
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of rat kidney...



Immunofluorescence analysis of mouse kidney...



Immunofluorescence analysis of rat kidney...