



Product Datasheet

Aquaporin 4 antibody (orb538467)


Description	Rabbit Polyclonal antibody to Aquaporin 4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Aquaporin-4 (AQP4)(NP_001641.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35KDa
Uniprot ID	P55087
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse brain tissue.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of human colon tissue.