

## Product Datasheet

### TRIM5 antibody (Biotin) (orb53850)

|                            |                                                                           |
|----------------------------|---------------------------------------------------------------------------|
| <b>Description</b>         | TRIM5 antibody (Biotin)                                                   |
| <b>Species/Host</b>        | Rabbit                                                                    |
| <b>Reactivity</b>          | Human                                                                     |
| <b>Conjugation</b>         | Biotin                                                                    |
| <b>Tested Applications</b> | ELISA                                                                     |
| <b>Immunogen</b>           | Recombinant Human Tripartite motif-containing protein 5 protein (1-493AA) |
| <b>Target</b>              | TRIM5                                                                     |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol            |
| <b>Form/Appearance</b>     | Liquid                                                                    |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.             |
| <b>Note</b>                | For research use only                                                     |
| <b>Isotype</b>             | IgG                                                                       |
| <b>Clonality</b>           | Polyclonal                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q9C035</a>                                                    |
| <b>Dilution Range</b>      | ELISA                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                           |