



Product Datasheet

ATP2A3 antibody (orb538545)


Description	Rabbit polyclonal antibody to ATP2A3
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ATP2A3(NP_005164.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	114kDa
Uniprot ID	Q93084
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human uterine tissue.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human brain tissue.