



## Product Datasheet

### BAG3 antibody (orb538574)


|                            |                                                        |
|----------------------------|--------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to BAG3                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                      |
| <b>Conjugation</b>         | Unconjugated                                           |
| <b>Tested Applications</b> | ICC, IF, IHC, IP, WB                                   |
| <b>Immunogen</b>           | Recombinant fusion protein of human BAG3(NP_004272.2). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.       |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.            |
| <b>Note</b>                | For research use only                                  |
| <b>Isotype</b>             | IgG                                                    |
| <b>Clonality</b>           | Polyclonal                                             |
| <b>MW</b>                  | 80kDa                                                  |
| <b>Uniprot ID</b>          | <b>O95817</b>                                          |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:100, IP: 1:50-1:100    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                        |



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human lung tissue.