


Product Datasheet

C6orf25 antibody (orb538668)

Description	Rabbit polyclonal antibody to C6orf25
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human C6orf25(NP_612121.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	O95866
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...