



Product Datasheet

CASC5 antibody (orb538714)


Description	Rabbit polyclonal antibody to CASC5
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human CASC5(NP_733468.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	300kDa
Uniprot ID	Q8NG31
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human CO...



Immunohistochemical staining of rat test...



Western blot analysis of extracts of var...