



Product Datasheet

CD34 antibody (orb538763)


Description	Rabbit polyclonal antibody to CD34
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human CD34(NP_001020280.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105KDa
Uniprot ID	P28906
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HUVEC cells.



Immunohistochemical staining of rat kidney tissue.



Immunohistochemical staining of mouse kidney tissue.