



## Product Datasheet

### CD36 antibody (orb538765)


|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CD36                    |
| <b>Reactivity</b>          | Human, Mouse, Rat                                     |
| <b>Conjugation</b>         | Unconjugated                                          |
| <b>Tested Applications</b> | ICC, IF, WB                                           |
| <b>Immunogen</b>           | A synthetic peptide of human CD36/SR-B3(NP_000063.2). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.      |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.           |
| <b>Note</b>                | For research use only                                 |
| <b>Isotype</b>             | IgG                                                   |
| <b>Clonality</b>           | Polyclonal                                            |
| <b>MW</b>                  | 70-100kDa                                             |
| <b>Uniprot ID</b>          | <b>P16671</b>                                         |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                       |



Immunofluorescence analysis of THP-1 cells...



Western blot analysis of extracts of various sources...



Western blot analysis of extracts of THP-1 cells...