

Product Datasheet

HBV pre S1/S2 protein antibody (orb5388)

Description	HBV pre S1/S2 protein antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from HBSAG_HBVD2 pre S1 protein (31-100/389 aa)
Target	Hepatitis B Virus Surface Antigen
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31/44 kDa
Uniprot ID	P03140
Expiration Date	12 months from date of receipt.