



Product Datasheet

CDK5R1 antibody (orb538810)

Description	Rabbit polyclonal antibody to CDK5R1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CDK5R1(NP_003876.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25KDa
Uniprot ID	Q15078
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SH-SY5Y C...



Western blot analysis of extracts of var...