



Product Datasheet

CDK5RAP1 antibody (orb538811)


Description	Rabbit polyclonal antibody to CDK5RAP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human CDK5RAP1(NP_057492.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68kDa
Uniprot ID	Q96SZ6
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of L-929 cell...



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...